



US006420181B1

(12) **United States Patent**  
**Novak**

(10) **Patent No.:** **US 6,420,181 B1**  
(45) **Date of Patent:** **\*Jul. 16, 2002**

(54) **FIELD MICROSPOT TEST METHOD FOR ON-SITE CHEMICAL TESTING**

(75) Inventor: **Thaddeus John Novak**, Bel Air, MD (US)

(73) Assignee: **The United States of America as represented by the Secretary of the Army**, Washington, DC (US)

(\*) Notice: Subject to any disclaimer, the term of this patent is extended or adjusted under 35 U.S.C. 154(b) by 0 days.

This patent is subject to a terminal disclaimer.

(21) Appl. No.: **09/296,602**

(22) Filed: **Apr. 23, 1999**

**Related U.S. Application Data**

(63) Continuation-in-part of application No. 08/763,181, filed on Dec. 11, 1996, now Pat. No. 5,935,862.

(51) Int. Cl.<sup>7</sup> ..... **G01N 33/00**; G01N 30/90

(52) U.S. Cl. .... **436/104**; 436/162; 436/164; 436/166; 436/169; 210/658

(58) Field of Search ..... 210/656, 658, 210/634, 635; 422/56, 58, 61, 68.1, 70; 436/103-105, 161-164, 166, 169, 808

(56) **References Cited**

**U.S. PATENT DOCUMENTS**

4,346,001 A \* 8/1982 Tyihak et al. .... 210/658

4,428,908 A \* 1/1984 Ashley et al. .... 422/61

\* cited by examiner

*Primary Examiner*—Lyle A. Alexander

(74) *Attorney, Agent, or Firm*—Ulysses John Biffoni; William W. Randolph

(57) **ABSTRACT**

A method, system and kit for detecting the presence of an analyte includes placing a solution containing the analyte in a microcapillary tube and placing the microcapillary tube in contact with a layer of sorbent material so that the solution is withdrawn from the microcapillary tube by capillary action. The sorbent material and solvent for the solution are selected so that the solvent is absorbed into the sorbent material and the analyte is adsorbed by the sorbent material and concentrated at the spot where the microcapillary tube contacts the sorbent material. A detector reagent is applied to the sorbent material to indicate the presence of the analyte.

**16 Claims, 2 Drawing Sheets**




FIG. 1







FIG. 3

FIG. 2



*FIG. 4**FIG. 5**FIG. 6*

1

## FIELD MICROSPOT TEST METHOD FOR ON-SITE CHEMICAL TESTING

This is a Continuation-In-Part of U.S. application Ser. No. 08/763,181, filed Dec. 11, 1996 now U.S. Pat. No. 5,935,862.

### GOVERNMENT INTEREST

The invention described herein may be manufactured, licensed, and used by or for the U.S. Government.

### BACKGROUND OF THE INVENTION

#### 1. Field of the Invention

This invention relates to a field test kit and method of on-site testing for the presence of contaminants and chemicals, and more particularly to a micro spot method for detecting the presence of a variety of chemicals and environmental contaminants.

#### 2. Description of the Prior Art

In view of biological hazards associated with toxic chemicals and environmental contaminants, regulations have been established by legislatures and environmental agencies to monitor a wide variety of chemicals and their byproducts. As a result, it is often necessary to conduct on-site inspections and analyses of various chemical spills, dump sites, and manufacturing facilities to detect environmental contaminants, hazardous conditions and to assure compliance with environmental regulations.

Advantages of on-site inspection and analysis of chemical sites include resolving ambiguities during the inspection, reducing the potential for contamination and cross-contamination of samples during travel to off-site testing laboratories, and providing a convenient method of performing a large number of preliminary tests to detect and screen for chemical contaminants. On-site inspection also provides a rapid indication of those samples which may possibly contain compounds that must be identified using more sophisticated laboratory analytical techniques. On-site testing also allows the level of concentration and spread of contamination from chemical spills to be readily surmised.

Reagent-based chemical detection and chromatographic methodologies are attractive for on-site testing and screening because many tests can be run in a short period of time and they are capable of providing visual presumptive evidence of the presence of a chemical substance in a sample. One methodology comprises classical spot tests that are normally carried out in depressions or wells of a porcelain spot plate. Conventionally, small amounts of a solution, which may contain chemical contaminants, are placed in the wells of the spot plate. Small quantities of different reagents are then added to the solution samples and a positive test is normally signified by a color change in the well of the spot plate. An advantage with these tests is that a number of tests can be carried out on a single plate. For example, as many as 12 different spot tests can be carried out on a small 3.5×4.5 inch spot plate. Another advantage is that it is possible to rapidly screen a large number of samples during a short period of time. However, as the concentration of the chemical substances become more dilute, it becomes more difficult to reliably detect the presence of the chemical substances. In most cases, the lower limit of detection is in the 1–100 microgram range.

Another methodology for screening samples and detecting target analytes in samples is use of thin layer chromatography (or TLC), which conventionally utilizes a plate

2

having a surface layer formed of a sorbent material or gel. In order to separate the components of the analyte obtained from a sample, a drop of solution is carefully applied above the bottom edge of a thin layer chromatography plate. Solutions suspected of containing target analytes are preferably deposited onto the surface of a TLC plate in the form of a drop to avoid a streaking pattern that would result if the device for applying the drop actually contacts the surface of the plate while a sample solution is being deposited. After the solvent evaporates, the residue on the plate is eluted with another solvent or solvent mixture (also known as the eluant) thereby causing the chemical components of the sample to migrate towards the top or opposite edge of the plate. When the proper conditions and eluant are chosen, each analyte migrates across the plate at a rate that is different from the other analytes. The elution step results in the different analytes separating from each other and settling at different regions and as diffuse spots along the path of migration. After the elution step, the plates are allowed to dry and then they are sprayed with a solution of visualization reagent (detector reagent). A persistent concern with thin layer chromatography is that the elution step of waiting for the solvent to completely wet the plate and for the analytes to migrate and separate is relatively time-consuming. In many instances, proper completion of the elution phase may exceed an hour and warrant involved techniques and quality control steps to assure adequate separation of the different analytes. Another concern involves situations where the analytes are present in such low concentrations that the detection signals obtained in the tests are weak and can possibly be misread. In summary, with thin layer chromatography the analytes in the sample migrate and separate into localized regions, as opposed to concentrated at spots or points.

### SUMMARY OF THE INVENTION

The micro spot test system and methodology of the present invention relates to an apparatus and method for the on-site testing of analytes contained in a sample by dissolving the analytes in a solvent and utilizing capillary deposition techniques to concentrate the analytes on sorbent materials. Detection sensitivity and accuracy for a range of concentrations of analytes is provided by applying a solution containing the analytes to a sorbent layer by capillary deposition so that the analytes in the solution become concentrated at the particular spot or point of deposition on the sorbent layer. The solutions are deposited by placing small diameter tubes containing the analyte solution in contact with the surface of the sorbent material so that the solutions are drawn from the small diameter tubes by capillary action. A detector reagent is then utilized to detect the presence of the analytes that are concentrated at the spot where the small diameter tube contacts the sorbent layer.

A system for chromogenically detecting the presence of chemical analytes includes a means for obtaining a sample solution containing the analytes; a device for the capillary deposition of the sample solution; chromatographic sorbent materials; and chromogenic detector reagents. Storage devices may be provided for the samples and for sample solutions; capillary deposition devices; the chromatographic sorbent materials and the chromogenic detector reagents.

Accordingly, one object of the present invention is to provide a compact chemical screening apparatus which is of a self-contained, efficient design for rapid screening of solutions for the presence or absence of target analytes.

Another object of the present invention is to provide a chemical screening device which is relatively simple to use

3

for sample solutions containing a wide range of analytes in a wide range of concentration levels.

These, together with still other objects of the invention, along with the various features which characterize the invention, are pointed out[]with particularity in the appended claims.

#### BRIEF DESCRIPTION OF THE INVENTION

Other objects and advantages of the invention will become apparent upon reading the following detailed description with reference to the attached drawings, wherein:

FIG. 1 is a plan view of a field test kit for performing on-site chemical analysis;

FIG. 2 is a plan view of a pipet with a micro-tip;

FIG. 3 is a plan view of removable micropipet tips;

FIG. 4 is a view of the end portion of a capillary tube;

FIG. 5 is a view generally showing the end portion of a capillary tube in contact with a sorbent layer; and

FIG. 6 is a plan view of a sorbent layer depicting the results of a micro spot test.

#### DETAILED DESCRIPTION OF THE PREFERRED EMBODIMENT

The micro spot system and methodology of detecting the presence of target analytes in a sample comprises the application of a solution containing the analytes to a chromatographic sorbent material by capillary action and adding a sufficient amount of a chromogenic detector reagent to form a chromogenic indicator when a target analyte is present in the sample. Apparatus for accomplishing this is generally shown in FIG. 1. In FIG. 1, the apparatus or kit 11 includes a bag or container 12 for storing the components of the system; at least one thin layer chromatographic plate or TLC plate 13; collecting devices such as cloth wipes 18 or swabs 19 for wiping surfaces for chemical residues; solvent containers 22; containers 23 for receiving the swabs and solvent solutions; reagent containers 24; and small diameter capillary or microcapillary tubes 25. A system or kit 11 for on-site detection and screening of a broad array of both volatile and non-volatile chemicals generally contains a wide variety of chemical reagents preferably stored in dry condition and in an inert atmosphere in small 2 to 3 ml. sealed containers or bottles 24, as represented in FIG. 1. Solvents for the samples and the reagents are stored in separate bottles 22 and an appropriate solvent is added to dissolve a reagent in the bottle in which it had been stored or in a separate container shortly prior to use. To insure a long shelf-life stability, the solid state reagents are preferably stored in dry condition, in an atmosphere that is free of moisture, and in which air has been displaced by an inert gas such as nitrogen or argon.

For purposes of this application, the term sample is defined as a representative fraction of the material that is to be processed and tested to detect the presence of an analyte. The sample may be a solid, such as soil, a liquid, such as water taken from a lake, or a vapor, such as fumes obtained from a chemical plant. An analyte is a chemical substance present in the samples that are being tested or analyzed. A solution is a homogeneous liquid that contains dissolved chemical substances. The analyte is a solute, which is defined as a chemical substance or mixture of chemical substances that dissolves in a solvent or a mixture of solvents to form a solution. A sample solution is a homogeneous liquid that contains dissolved chemical substances (i.e., the analytes or solutes) and which is derived by

4

washing, extracting, or eluting a sample with a solvent or mixture of solvents. For example, surface wipes 18 or swabs 19 of polyester or similar material are used to obtain a sample by wiping a suspected surface. A sample solution is obtained for analysis by washing, extracting or eluting the wipe in a container 23 with a suitable solvent such as acetone, dichloromethane, hexane, etc. Soil samples can be washed, extracted or eluted in separate containers to obtain sample solutions. Aqueous samples suspected of containing a target analyte can be extracted with an immiscible solvent which is capable of extracting the analytes believed to be therein. In addition, solid phase extraction (SPE) or solid phase microextraction (SPME) techniques can be used to extract analytes from water for analysis using the micro spot tests.

Once the solution or liquid extract has been formed, and where necessary the extract has been concentrated by evaporation, a tube with a small diameter bore or opening 25, such as a small diameter capillary or microcapillary tube is used to collect and dispense small amounts of the solution onto the surface of plate 13 by capillary action. Preferably, the plates 13 are thin layer chromatographic plates or TLC plates having a surface layer formed of a chromatographic sorbent material. A sorbent material is a material that has both absorption and adsorption characteristics. Absorption is defined as the penetration of liquids into the bulk of a porous material somewhat like a sponge soaking up water. Adsorption is a process whereby a chemical substance, an analyte, sticks, clings or adheres to the surface of a solid constituent, the adsorbent. In FIG. 1, the plates have been provided with scoring lines 14 to divide the plate 13 into a plurality of separate sections 15 that serve as different test sites. Generally, the amount of sample delivered to a test site on the chromatographic material from a microcapillary tube having a length of one and one quarter inches is on the order of from about 0.1 microliters (for an approximate 0.05 mm diameter microcapillary opening) to about 30 microliters (for an approximate 1.6 mm diameter microcapillary opening) of sample. In most instances sample size will be on the order of from about 0.5 microliters (for an approximate 0.1 mm diameter microcapillary opening) to about 5 microliters (for an approximate 0.4 mm diameter microcapillary opening) and preferably, the sample will be on the order of from about 1 microliter (for an estimated 0.2 mm diameter microcapillary opening) to about 3 microliters (for an estimated 0.25 mm diameter microcap opening). Microcapillary tubes having longer lengths can be used. If desired, the microcapillary tube 25 can be held with commercially available holders or forceps.

The term "microcapillary tube" includes any tube made from glass, plastic or other material having a small diameter opening that is capable of dispensing liquid from (or drawing liquid into) the opening by capillary action. Examples of small diameter capillary tubes are those marketed by Drummond and sold under the trademarked name of Microcaps. Another type of tube having a small diameter opening is a micropipet. A micropipet is a glass or plastic tube having a small diameter opening (or capillary opening) at one end and an enlarged opening at the other end of the micropipet, as generally shown by micropipet 27 in FIG. 2. Examples of micropipets are Micro-tip polyethylene pipets sold by Micro Mole Scientific. One benefit of a micropipet is that if the top of the bulb is cut off, as shown in FIG. 2, the larger end functions as a funnel for holding a larger volume of fluid sample than could normally be held or drawn into a capillary tube. Consequently, a larger volume of sample (such as 10 microliters or more) can be used to achieve a higher detec-

5

tion sensitivity with respect to the concentration of analyte that can be detected. Replaceable micropipet tips **31**, as shown in FIG. **3**, are examples of additional devices that have small diameter openings. A removable micropipet tip **31** would be placed on the end portion of another tube **33** or container so that the liquid in the tube or container would be withdrawn by capillary action when placed in contact with a chromatographic sorbent material. An additional deposition control could be achieved by use of an in-line filter element **32** in a microcapillary device. One example of a micropipet tip is available under the tradename of Plasti-brand autoclavable nonsealing filter tips and another example of a micropipet tip without a filter is Catalog No. 71-6311-10 from PGC Scientific, Gaithersburg, Md. An example of a micropipet tip with a filter is Catalog No. 71-6311-16 from PGC Scientific, Gaithersburg, Md. While microcapillary tubes, micropipets and micropipet tips have been distributed for use with bulbs or other devices for forcing liquid out of the tubes, use of such pressure devices for forcing liquid from microcapillary tubes is contrary to the methodology of the present invention where the solution containing the analyte is deposited by capillary action. Further, while some methodology for applying a sample to a TLC plate with a microcapillary tube includes moving the tip of the microcapillary tube as the sample is being applied, use of such methodology that includes moving the tip of the microcapillary tube is contrary to the methodology of the present invention where all of the analyte contained in the sample solution must be adsorbed in the smallest volume of sorbent.

To avoid breaking extremely thin microcapillary tubes, micropipets and micropipet tips during use as the ends portions of these devices contact sorbent surfaces, it is possible to use various holding devices such as forceps and small clamps. While microcapillary tubes have been found acceptable for most applications, where conditions or technique warrant, the end portions of the tubes can be formed with thickened wall portions as shown in FIG. **4**, where the thickness of the wall portions **35** are at least equal to the diameter of the opening **36** of the tube **25**. Increasing the wall thickness to at least twice the diameter of the opening not only strengthens the end portion of the tubes for adverse use conditions, but also provides a larger contact surface area relative to the size of the opening and thereby promotes higher circumferential contact and seal of the microcapillary tubes with the sorbent material.

In general, the concentration of analyte in solution that is capable of being detected using a microcapillary tube to apply the solution containing the analyte to a thin layer chromatography plate (or TLC plate) in a micro-spot test is inversely related to the total volume of solution applied to the plate. For example, analytes present in solutions at low concentration levels can be detected by increasing the total volume of solution applied by capillary action to a TLC plate in the micro spot test. When using a microcapillary tube such as a micropipet to apply a relatively large volume of solution to a TLC plate, the sample solution being drawn from the tip of the micropipet should initially come in contact only with the area directly beneath the opening of the micropipet that contacts the TLC plate. One way for this condition to be met is that the solution can be added to the large end of the micropipet tube in small aliquots, so that the solution wets the TLC plate by capillary action. Another way for this condition to be met is that when a large volume of solution is added in a single aliquot to the micropipet, sufficient pressure should be exerted so that the micropipet tip completely contacts the TLC plate. When this occurs, the solu-

6

tion will wet the plate by capillary action. To promote well-defined small spots, the liquid solution should not be allowed to leak or flow from the juncture of the micropipet tip with the surface of the chromatographic sorbent layer of the TLC plate. A further methodology is to use a micropipet tip that contains a filter or other device within the tip that slows the flow of solution to ensure that the solution wets the TLC plate by capillary action when the solution is applied to the TLC plate by placing the tip in contact with the sorbent layer of the TLC plate.

Preferably, the sorbent layer should be capable of acting both as an adsorbent and as an absorbent. It is believed that when the sample solution, which contains an analyte (i.e., the solute) dissolved in a solvent, is applied to the sorbent by capillary action using a microcapillary device, the solute separates from the solution because it adheres or clings to the walls of the pores in a small volume of the sorbent immediately surrounding the point of application of the sample solution due to adsorption, while the solvent, a fluid which consequently has been freed from the solute, fills the voids in the pores of the sorbent due to absorption. This phenomenon results in the analyte concentrating within the sorbent layer and being localized in a small volume of the sorbent or "spot", while the solvent freely wets a substantial volume of sorbent. Since the volume of sorbent in which the solute is adsorbed is a small fraction of the volume in which the solvent is absorbed, the analyte becomes highly concentrated and consequently, high sensitivity of detection is made possible in the micro spot tests. When the analyte is present in very low concentration levels, application of the sample solution to the sorbent layer will result in the analyte concentrating in a very small volume of sorbent, and hence, will produce only a very small spot. Whereas, if the analyte is present in a somewhat higher concentration level, application of the sample solution to the sorbent layer will result in the analyte concentrating in a somewhat larger volume of sorbent, and hence, will produce a somewhat larger spot. FIG. **5** depicts a view where the end portion of a microcapillary tube has been placed in sufficient contact with a sorbent layer so that as the solution containing the analyte leaves the opening in the end portion of the tube, the analyte is adsorbed in a small localized region or spot about the point where the tip or end of the microcapillary tube contacts the sorbent layer, while the solvent spreads throughout the porous medium as it wets and is absorbed into the sorbent layer.

The chromatographic sorbent material is preferably a thin-layer chromatography (or TLC) plate. TLC plates are commonly found containing a silica gel or alumina coating. One example is MK6F Silica Gel 60A TLC plates, Catalog No.4861-110 from Whatman, Inc., Clifton N.J. 07014, which contains a 250 micron thick layer on a 1 by 3 inch glass microscope slide is a preferred solid support for the micro spot tests. These plates are relatively easy to handle and they contain a solid substance that fluoresces brightly when illuminated with short wavelength UV light. Therefore it is possible to detect substances that absorb UV light by "fluorescence quenching", as well as by other detection and visualization methods. A non-limiting list of suitable TLC strips which can be used in carrying out the invention include Diamomd MK6F Silica Gel 60A TLC plates, Catalog No. 4500-100 from Whatman, Inc., Clifton N.J. 07014, which contains a 250 micron thick layer on a 1 by 3 inch glass microscope slide; Silica Gel HL, 250 micron thick layer, Cat. No. 46931, Analtech, Newark, DE 19714; Silica Gel HLF, 250 micron layer, Cat. No. 47931, Analtech, Newark, DE 19714; Ammonium Sulfate (5%) Modified

Silica Gel H, 250 micron layer, Self Charring Plates, Cat. No. 74031 (without indicator) and Cat. No. 75031 (with fluorescent indicator); Silica Gel F-254 TLC media, plastic backed, layer thickness 0.25 mm, Cat # 5775 from E. M. Laboratories, Elmsford, N.Y. 10523; Silica Gel F-254 TLC media, aluminum backed, layer thickness 0.2 mm, Cat # 5539 from Alltech Associates, Deerfield, Ill. 60115; Silica Gel TLC media, plastic backed, layer thickness 100 microns, Product Number 13179, Cat. # 4G 6801, Eastman Kodak Co., Rochester, N.Y. 14650; C<sub>18</sub>/Silica Gel, 250 micron thick layer, Cat. No. 17021, Analtech, Newark, Del. 19714; NH<sub>2</sub>/Silica Gel, 250 micron thick layer, Cat. No. 18021, Analtech, Newark, Del. 19714; CN/Silica Gel, 250 micro thick layer, Cat. No. 19021, Analtech, Newark, Del. 19714; Nano-SIL G High Performance Thin-Layer Chromatography (HPLTC) Plates, Catalog 81841, Alltech, Inc., Deerfield, Ill. 60015; Nano-SIL-NH<sub>2</sub>/UV (Amino). Catalog No. 8100026, Alltech, Inc., Deerfield, Ill. 60015; Nano-SIL-CN/UV (Cyano) Catalog No. 8110022, Alltech, Inc., Deerfield, Ill. 60015; Reversed Phase Sil Gel 60, RP-2 (Dimethyl bonded) Cat. No. 5746, RP-8 (Octyl bonded) Cat. No. 15388-7, RP-18 (Octadecyl bonded) Cat. No. 15389-7, Alltech Inc., Deerfield, Ill. 60015; "hybrid plates" (one plate designed for both reverse-phase and normal phase TLC), Catalog Number 818144, Alltech, Inc., Deerfield, Ill. 60015; Avicel Microcrystalline Cellulose Uniplates, 250 micron thick layer, Cat. No. 05061 (without indicator) and Cat. No. 06061 (with fluorescent indicator), Analtech, Newark, Del. 19714; SILCEL-Mix 25 UV254, Catalog No. 810043, Alltech, Inc., Deerfield, Ill. 60015; ALOX-100 UV254, Catalog No. 807033, Alltech, Inc., Deerfield, Ill. 60015; GUR N-25 UV254, Catalog No. 810073, Alltech, Inc., Deerfield, Ill. 60015; Nano-SIL C18-100 UV254, Catalog No. 811062, Alltech, Inc., Deerfield, Ill. 60015; SIL N-HR/UV254, Catalog No. 804023, Alltech, Inc., Deerfield, Ill. 60015; CEL 300 AC-30%, Catalog No. 801043, Alltech, Inc., Deerfield, Ill. 60015; CEL 300 DEAE, Catalog No. 801073, Alltech, Inc., Deerfield, Ill. 60015; Polyamide 6 UV254, Catalog No. 803023, Alltech, Inc., Deerfield, Ill. 60015; ALOX N/UV254, Catalog No. 802021, Alltech, Inc., Deerfield, Ill. 60015; Instant Thin Layer Chromatography Polysilicic Acid Gel Impregnated Glass Fiber Sheets with Fluorescent Indicator, Product Number 51435, Gelman Instruments, Ann Arbor, Mich. 48106; Instant Thin Layer Chromatography Sheets, Type SG, Product Number, 61886, Gelman Instrument Co., Ann Arbor, Mich. 48106; TLC Plates, Silica Gel 60 F-254, aluminum backed, layer thickness 0.2 mm, Product # 37360, Catalog # Z19,329-1, Aldrich Chemical Co., Milwaukee, Wis. 53233; Silica Gel IB Flexible (plastic backed) Sheets for Thin Layer Chromatography, layer thickness 250 microns, Product Number 4462-02, J. T. Baker, Inc., Phillipsburg, N.J. 08865; Aluminum Oxide IB Flexible (plastic backed) Sheets for Thin Layer Chromatography, layer thickness 200 microns, Product Number 4466-00, J. T. Baker, Inc., Phillipsburg, N.J. 08865; Reversed Phase (hydrocarbon impregnated) HPTLC Uniplates, 150 micron thick layer, Cat No. 54377 (without indicator) and Cat. No. 55377 (with fluorescent indicator), Analtech, Newark, Del. 19714; MKC18F Reversed Phase TLC plates, glass backed (1"x3" plates), layer thickness 200 microns, Cat. # 4803-110 from Whatman, Inc., Clifton, N.J. 07014; H-RP2F (ethyl bonded silica gel) Reversed Phase TLC plates, layer thickness 50 microns, Cat. No. 08527, Analtech, Newark, Del. 19714; Polyram Ionex 25 SA-NA Ion Exchange Resin and Silica Gel Mixed Layer on Plastic, Catalog Number M806013, Bodman Chemical Co., Aston, Pa. 19014; Polygram Ionex 25 SB-AC Ion Exchange Resin

and Silica Gel Mixed Layer on Plastic, Catalog Number 806023, Bodman Chemical Co., Aston, Pa. 19014; and 2000 micron thickness Silica Gel G Preparative Uniplates, Catalog Number 01055, Analtech, Inc., Newark, Del. 19714. The composition of adsorbent coatings contained on the listed TLC plates include silica gel, high performance thin layer chromatography (HPTLC) silica gel, polysilicic acid, aluminum oxide, cellulose, polyamide, reversed phase silica Gel C<sub>2</sub> (dimethyl bonded), reversed phase silica gel C<sub>2</sub> (ethyl bonded), reversed phase silica gel C<sub>8</sub> (octyl bonded), reversed phase silica gel C<sub>18</sub> (octadecyl bonded), acetylated cellulose, silica gel modified with amino groups, silica gel modified with cyano groups, Kieselghur impregnated with hydrocarbons, anionic and cationic anion exchange resins, diethylaminoethyl cellulose, and mixtures of the listed sorbents.

If the thin-layer chromatography (or TLC) sheets or plates are those which are commercially available, they can be further scored into small sample areas within the sheet, i.e. 1/2 inch x 1/2 inch or a similar size. The scoring of the plate reduces the likelihood that the liquid detector reagent applied in one spot test will creep into the sections reserved for other spot tests. Alternatively, the method can be carried out using TLC sheets that are specifically made to carry out the micro spot tests of the present invention.

After the sample solution is deposited on a TLC plate, a short period of time, generally on the order of about one minute, is allowed to elapse during which the solvent evaporates. Then, a sufficient amount of a detector reagent is added to the location on the TLC plate where the solution containing the analyte was deposited. It is to be understood that the solution containing the analyte can be deposited on a plurality of test sites on the TLC plates and different detector reagents can be used to test for the presence of target analytes or functional groups. Alternately, different sample solutions can be applied each to a different spot on a single TLC plate, and one detector reagent can be added to all of the analyte spots on the plate.

Examples of chromogenic detector reagents include bromocresol green; 7,7,8,8-tetracyanoquinodimethane (TCNQ); gold chloride (without NaOH); gold chloride/NaOH solution (i.e., the slash between the reagents means "followed by" and it is applicable to all of the reagent combinations); 4-(4'-nitrobenzyl)pyridine/NaOH; cholinesterase/indoxyl acetate; sodium pyrophosphate peroxide/aromatic amine; potassium bismuth iodide; 1,3-diisonitrosoacetone guanidinium salt; bis(diethylamino)benzophenone oxime; bis(diethylamino)benzophenone; bis(dimethylamino)thiobenzophenone; phenylazoformic acid 2-diphenylhydrazide; diphenylcarbazon; diphenylthiocarbazon; mercuric salt; diethyldithiocarbamic acid silver salt; 2,2'-dithiobis(5-nitropyridine); 5,5'-dithiobis(2-nitrobenzoic acid) i.e., Ellman's Reagent; molybdenum oxide in sulfuric acid; ammonium molybdate; iodine/starch; and sulfuric acid (4M). Below is a non-limiting partial list of other detector reagents that can be used in the methodology described herein. Other detector reagents are contained in references 1-38 which are incorporated herein by reference.

TABLE 1

Partial List of Detector Reagents for Microspot Tests	
Detector Reagent	Target Analytes
alizarin	cations
aluminum chloride	flavonoids

TABLE 1-continued

Partial List of Detector Reagents for Microspot Tests	
4-aminoantipyrine/potassium hexacyanoferrate	phenols
o-aminodiphenyl/phosphoric acid	sugars
4-aminohippuric acid	reducing sugars
o-aminophenol/phosphoric acid	sugars
ammonia	tetracyclines
ammonium cerium(IV)nitrate	polyalcohols
ammonium cerium(IV)nitrate/nitric acid	alpha-hydroxy acids,
	alpha-keto acids,
	mercaptans
ammonium cerium(IV)sulfate	alkaloids
ammonium iron(III)sulfate	flavonoids
ammonium iron(III)sulfate	alkaloids
ammonium molybdate/crystal violet	phosphoric acid
ammonium molybdate/tin(II)chloride	phosphoric acids
ammonium thiocyanate/iron(II)sulfate	peroxides
aniline/diphenylamine/phosphoric acid	reducing sugars
aniline/phosphoric acid	sugars
aniline phthalate	reducing sugars, anions
	of halogen oxyacids
anisaldehyde/sulfuric acid	sugars, steroids, terpenes
p-anisaldehyde	reducing sugars
p-anisidine phthalate	reducing sugars
anthrone	ketoses
antimony(III)chloride	flavonoids
antimony(III)chloride	vitamin A and D,
	carotenoids, steroids,
	sapogenins, steroid
	glycosides, terpenes,
	sapogenins resins,
	steroid
aurin tricarboxylic acid (aluminon)	aluminum ions,
	chromium ions,
	lithium ions
2,2'-bipyridine/iron(II)chloride	phenols, vitamin E,
	reducing compounds
bismuth chloride	sterols
boric acid/citric acid	quinolines
bromine/fluorescein/silver nitrate	insecticides
bromocresol green (or bromocresol green)	organic and inorganic
	acids,
	organic acids
bromocresol green/bromophenol blue/	dicarboxylic acids,
potassium permanganate	halogen ions
bromocresol purple	phenols
bromophenol blue/methyl red	lipids, sulfur compounds
bromosuccinimide/fluorescein	vitamin C
cacotheline	sugars
carbazole/sulfuric acid	polysaccharides
carmine	organic and inorganic
cerium(IV)sulfate	iodine compounds
cerium(IV)sulfate/nitric acid	polyphenylenes
cerium(IV)sulfate/sulfuric acid	alkaloids, iodo-organic
	compounds
chloramine-T	caffeine
chloramine-T/trichloroacetic acid	digitalis glycosides
1-chloro-2,4-dinitrobenzene	nicotinic acid,
	nicotinamide, pyridoxol
chlorosulfonic acid/glacial acetic acid	triterpenes, sterols,
	steroids
chromosulfuric acid	organic compounds
chromotropic acid	methylenedioxyphenyl-
	type compounds nicotine,
	hydrastine, sesamine
cinnamaldehyde/acetic anhydride/sulfuric acid	steroid sapogenins
cinnamaldehyde/hydrochloric acid	indole derivatives
cobalt(II)chloride	organic phosphate esters
cobalt(II)/lead nitrite	ammonium ions,
	potassium ions
cobalt(II)nitrate/ammonia	barbiturates
cobalt(II)nitrate/lithium hydroxide	barbiturates
cobalt(II)thiocyanate	alkaloids, amines
copper acetate/potassium hexacyanoferrate(II)	higher fatty acids
copper acetate/rubeanic acid	higher fatty acids
copper chloride	oximes
copper sulfate/quinine/pyridine	barbiturates,
	thiobarbiturates

TABLE 1-continued

Partial List of Detector Reagents for Microspot Tests	
5 copper(II)sulfate/sodium citrate	flavonoids,
	coumarins with
	o-dihydroxy groups
	straight-chain lipids
	desoxyribonucleosides
	2-deoxy-sugars
	aldehydes, ketones
	phenols
	saturated and
	unsaturated lipids
	free fatty acids
	2,7'-dichlorofluorescein/aluminum chloride/
	iron(III)chloride
10 alpha-cyclodextrin	
cysteine/sulfuric acid	
3,5-diaminobenzoic acid/phosphoric acid	
o-dianisidine	
2,6-dibromoquinone chlorimide	
2',7'-dichlorofluorescein (fluorogenic indicator)	
	2,6-dichlorophenolindophenol/silver nitrate
15 2,6-dichlorophenolindophenol sodium salt	
	2,6-dichloroquinone chlorimide
20 dicobalt octacarbonyl	
diethylamine/copper(II)sulfate	
diethyl malonate	
	dimedone/phosphoric acid
	4-dimethylaminobenzaldehyde/acetic acid/
25 phosphoric acid	
4-dimethylaminobenzaldehyde/acetylacetone	
4-dimethylaminobenzaldehyde/hydrochloric	
acid	
4-dimethylaminobenzaldehyde/sulfuric acid	
dimethylaminobenzylidenetherhodanine	
30 4-dimethylaminocinnamaldehyde	
dimethyl-p-phenylenediamine dihydrochloride	
dimethyl-p-phenylenediamine dihydrochloride/	
trichloroacetic acid	
1,3-dinitrobenzene	
35 3,5-dinitrobenzoic acid	
3,5-dinitrobenzoic acid	
2,4-dinitrofluorobenzene	
2,4-dinitrophenylhydrazine	
	3,5-dinitrosalicylic acid
40 diphenylamine	
diphenylamine/palladium(II)chloride	
diphenylamine/zinc chloride	
diphenylboric acid-beta-aminoethyl ester	
	diphenylcarbazide
45	
	diphenylcarbazone
50	
diphenylcarbazone	
diphenylpicrylhydrazyl	
2,5-diphenyl-3-(4-styrylphenyl)-tetrazolium	
chloride	
dipicrylamine	
55 dithizone	
4,4'-dithiodianils	
Dragendorff reagent	
60 ethylenediamine	
ethylenediamine/potassium hexacyanoferrate	
65 Fast blue salt B	
fluorescein	



TABLE 1-continued

Partial List of Detector Reagents for Microspot Tests	
fluorescein/ammonia	purines, pyrimidines, barbiturates
fluorescein/bromine	unsaturated compounds
fluorescein/hydrogen peroxide	hypnotics containing bromine
fluorescein/rhodamine-B/sodium carbonate	chlorinated hydrocarbons, heterocyclic compounds
formaldehyde/hydrochloric acid	indoles, indole derivatives
formaldehyde/phosphoric acid	steroid alkaloids, steroid sapogenins, phenothiazine derivatives
formaldehyde/sulfuric acid	aromatic compounds
furfural/sulfuric acid	carbamate esters
glucose/aniline	acids
glucose/phosphoric acid	aromatic amines
glyoxalbis(2-hydroxyanil)	cations
hydrazine sulfate	piperonal, vanillin, ethyl vanillin
hydrochloric acid	glycols
hydrogen peroxide	aromatic acids
4-hydroxybenzaldehyde/sulfuric acid	sapogenins, corticosteroids
hydroxylamine/iron(III)chloride	lactones, esters, amides, anhydrides of carboxylic acids
8-hydroxyquinoline	barium ions, strontium ions, calcium ions
8-hydroxyquinoline/hypobromite	guanidine derivatives
indandione	carotenoid aldehydes
iodine	general detection reagent
iodine/potassium iodide	alkaloids
iodine/potassium iodide	organic compounds
iodine/sulfanilic acid/N-(1-naphthyl)ethylene diamine	hydroxylamines
iodine/sulfuric acid	organic compounds containing nitrogen, polyethylene glycols, polyethylene glycol derivatives
iron(III)chloride	phenols, hydroxamic acids
iron(III)chloride/iodine	xanthine derivatives
iron(III)chloride/potassium hexacyanoferrate (III)/arsenite	thyroid hormones, iodine containing compounds
iron(III)chloride/sulfosalicylic acid	thiophosphate esters
iron(III)chloride/sulfuric acid	indol derivatives
iron(II)thiocyanate	peroxides
isatin/sulfuric acid	thiophene derivatives
isatin/zinc acetate	amino acids
isonicotinic acid hydrazide	ketosteroids
lead acetate (basic)	flavonoids
lead(IV)acetate	1,2-diol groups
lead(IV)acetate/rosaniline	1,2-diol groups
magnesium acetate	anthraquinone glycosides
mercury chloride/diphenylcarbazone	barbiturates
mercury(II)chloride/potassium iodide	steroid alkaloids
mercury(I)nitrate	barbiturates
methylene blue	sulfate esters of steroids
methylumbelliferone (fluorogenic detector reagent)	heterocyclic compounds
methyl yellow	containing nitrogen
molybdato-phosphoric acid	chlorinated insecticides
	reducing compounds, lipids, sterols, steroids
	aluminum ions
morin	sugars
1,3-naphthalenediol/phosphoric acid	sugars
1,3-naphthalenediol/sulfuric acid	sugars, uronic acids
1,3-naphthalenediol/trichloroacetic acid	guanidine derivatives
1-naphthol/hypobromite	amino acids,
naphthoquinone-sulfonic acid sodium salt	aromatic amines

TABLE 1-continued

Partial List of Detector Reagents for Microspot Tests	
1-naphthylamine	3,5-dinitrobenzoic acid esters, dinitrobenzamides
ninhydrin	amino acids, amines, amino-sugars, amines, amino acids
ninhydrin/cadmium acetate	amino acids
ninhydrin/copper(II)nitrate	amines
ninhydrin/tin(II)chloride	iron ions
2-nitroso-1-naphthol-4-sulfonic acid	thiophosphate esters, organophosphorus insecticides
palladium(II)chloride	sugars
phenol/sulfuric acid	reducing sugars
m-phenylenediamine	conjugated 3-ketosteroids
p-phenylenediamine/phthalic acid	dehydroascorbic acid
o-phenylenediamine/sulfuric acid	alpha-keto acids
o-phenylenediamine/trichloroacetic acid	germanium
phenylfluorone	dehydroascorbic acid
phenylhydrazine	sterols, steroids
phosphoric acid	digitalis glycosides
phosphoric acid/bromine	iron(III) ions
potassium hexacyanoferrate(II)	adrenaline and derivatives
potassium hexacyanoferrate(III)	vitamin B1
potassium hexacyanoferrate(III)/iron(III)chloride	reducing compounds, phenols, amines, thiosulfates, isothiocyanates
potassium hexacyanoferrate(III)/phosphate buffer	adrenaline
potassium hexacyanoferrate(III)/potassium hexacyanoferrate(II)	morphine
methanolic potassium hydroxide	coumarins, anthraquinone glycosides
potassium iodide/hydrogen sulfide	heavy metal ions
potassium iodide/starch	peroxides
potassium iodine plateate	alkaloids, other organic compounds
potassium permanganate (alkaline)	containing nitrogen, ketosteroids
potassium permanganate (neutral)	reducing compounds, sugars, aromatic polycarboxylic acids, polyalcohols
quinalizarin	easily oxidizable compounds
p-quinone	cations
resorcinol/zinc chloride/sulfuric acid	ethanolamine
resorcyaldehyde/sulfuric acid	phthalate esters
rhodamine 6G	16-dehydrosteroids
rhodanine	lipids
rhodizonic acid sodium salt	carotenoid aldehydes
rubeanic acid	barium ions, strontium ions
silver nitrate	lead ions, cobalt ions, copper ions
silver nitrate/ammonia	manganese ions, nickel ions, mercury ions, bismuth ions
silver nitrate/ammonia	phenols
silver nitrate/ammonia/fluorescein	sugars, sugar alcohols
silver nitrate/ammonia/sodium chloride	reducing substances
silver nitrate/ammonia/sodium methoxide	halogen ions
silver nitrate/bromophenol blue	thioacids
silver nitrate/fluorescein	sugars
silver nitrate/formaldehyde	purines
silver nitrate/potassium dichromate	alkylsulfonic acids
silver nitrate/potassium permanganate	arylsulfonic acids
silver nitrate/sodium dichromate	chlorinated insecticides, dieldrin, aldrin, lindane
silver nitrate/sodium hydroxide	barbiturates
	reducing compounds
	purines
	sugars, polyalcohols

TABLE 1-continued

Partial List of Detector Reagents for Microspot Tests	
sodium meta-periodate	hydroxyamino acids
sodium meta-periodate/4-nitroaniline	serine, threonine
sodium nitrite/hydrochloric acid	deoxy-sugars
sodium nitroprusside	indoles, thiazoles
	compounds with
	sulfhydryl group
sodium nitroprusside/acetalddehyde	secondary aliphatic
	and alicyclic amines
sodium nitroprusside/ammonia	hemlock alkaloids
sodium nitroprusside/hydroxylamine	thiourea derivatives
sodium nitroprusside/potassium hexacyano-	aliphatic nitrogen
ferrate(III)	compounds cyanamide,
	guanidine, urea,
	thiourea, thiourea
	derivatives, creatine,
	creatinine
sodium nitroprusside/potassium permanganate	sulfonamides
sodium nitroprusside/sodium hydroxide	methyl ketones
	activated methylene
	groups
sodium nitroprusside/sodium meta-periodate	deoxy-sugars
sodium pentacyanoamino ferrate(II)	urea, thiourea,
	guanidines
sodium sulfide (aqueous solution)	hydrogen sulfide group
sodium tetraphenylboron	alkaloids
sodium tetraphenylboron/rhodamine B	potassium ions
sodium thiosulfate/copper(II)acetate	antimony ions
starch	amylases
sulfanilic acid/1-naphthylamine	nitrosamines
sulfuric acid	general visualization
sulfuric acid/hypochlorite	digitalis glycosides
	reagent
tetracyanoethylene	aromatic hydrocarbons,
	phenols, heterocyclic
	compounds
tetranitrodiphenyl	cardiac glycosides
tetraphenyldiboroxide	flavones
tetrazolium blue	corticosteroids,
	reducing compounds
thiobarbituric acid	sorbic acid
thymol/sulfuric acid	sugars
thymol blue	dimethylamino acids
tin(II) chloride/hydrochloric acid/4-dimethyl-	aromatic compounds
aminobenzaldehyde	containing nitro groups
tin(II) chloride/potassium iodide	gold ions
tin(IV) chloride	triterpenes, sterols,
	steroids, phenols,
	polyphenols
titan yellow	cadmium ions
p-toluenesulfonic acid	steroids, flavonoids,
	catechins
o-tolidine (fluorogenic detector reagent)	chlorinated insecticides
toluidine blue	acidic polysaccharides
trichloroacetic acid	steroids, digitalis
	glycosides, Veratrum
	alkaloids, vitamin D
	steroids
trifluoroacetic acid	thiophosphate pesticides
N,2,6-trichloro-p-benzoquinoneimine	cardiac glycosides
2,4,6-trinitrobenzoic acid	reducing sugars,
2,3,5-tripheyltetrazolium chloride	corticosteroids,
	reducing compounds
tungstophosphoric acid	reducing compounds,
	lipids, sterols,
	steroids
urea/hydrochloric acid	sugars
vanillin/hydrochloric acid	catechins
vanillin/phosphoric acid	steroids
vanillin/potassium hydroxide	ornithine, lysine,
	proline, amines
vanillin/sulfuric acid	higher alcohols, phenols,
	steroids, essential oils
violuric acid	alkali and alkaline
	earth metal ions
xanthidrol	tryptophan,
	indole derivatives

TABLE 1-continued

Partial List of Detector Reagents for Microspot Tests	
5 zinc chloride	steroid sapogenins,
	steroids
zinc uranyl acetate	sodium ions
zirconyl chloride/alizarin/hydrochloric acid	fluorine ions
zirconyl chloride/citric acid	glycosides

10 It is generally expected that one drop of a detector reagent(s) will be sufficient in order to produce a result. However, in tests where two detector reagents are added in sequence, the second detector reagent should be added about 15 two minutes (or more) after the first reagent. In some tests heat can be used to accelerate the reaction of a detector reagent with the analyte. If the test is positive, in most cases a small spot within the reagent spot on the TLC sheet 20 changes color, often instantaneously, but with low analyte levels it will generally require a longer period of 15–60 seconds. Since the level of analyte tends to be directly related to the size of the color change within the spot, some quantification of the analyte levels may be possible. With the micro spot test, positive tests were obtained where the amount of analyte in a solution applied to a TLC plate was 25 as low as 10 ng (see examples 1–5 below). One analyte, EMPTA, was also detected at the 1 ng level (see example 3 below). It is also possible to detect an analyte where the detection signal produced with the disclosed methodologies is a change in the intensity of color rather than an actual 30 color change. For example, strong bases can be detected in disclosed micro spot test of the present invention where Bromcresol Green is the detector reagent. The detection signal that was produced was a small dark blue spot within 35 a light blue spot in the test with Bromcresol Green when the analyte was a strong base.

FIG. 6 is a plan view of a sorbent layer depicting how the results of a microspot test may appear after a detector reagent has been added to the surface of a sorbent layer. The two circular spots or regions 40 and 42 are generally representative of where the detector reagent solution has been added to a sorbent layer. The spot to the left 40 represents the results when a negative result is produced (no analyte is present) and the spot to the right 42 represents the results when a positive result is produced (an analyte is present, producing a chromogenic indication). The right spot is shown to contain a smaller spot or point 43 in the center of the detector reagent spot to represent that the unknown has been localized about the point where the microcapillary tube 50 contacts the sorbent material. While there may be some expansion of the spot relative to the size of the opening, such as would occur with a relatively high concentration of analyte or use of different solvents, the analyte generally remains concentrated at the identified spot of deposition.

55 It was found that high detection sensitivity is attained with the disclosed methodology where the analyte is dissolved in a solvent and then applied in small amounts to a thin-layer chromatography media by capillary action using a micro-capillary tube. Thus, the test method of the present invention is referred to as micro spot tests due to the minute quantities of analyte that are capable of being detected. It is further believed that the increased sensitivity of the micro spot tests is due at least in part to the fact that the analyte remains concentrated and localized at the spot of deposition.

65 The micro spot test method is generally solvent dependent, with respect to both the solvents for the analyte and the detector reagent. In general, the solvent for target

analytes should be selected so that the analytes are concentrated in a small spot when the solution containing the analytes is applied to a TLC plate with a microcapillary tube. The solvent for the detector reagent should be selected so that the spot containing the analytes is not enlarged too much or washed from the TLC plate when the detector reagent solution is applied. For example, the silica gel and aluminum oxide (alumina) sorbents in TLC plates are polar compositions. When a solution that contains an analyte is applied to a polar sorbent coating such as silica gel or alumina with a microcapillary tube, a smaller and more compact spot will tend to form if the solvent for dissolving the analyte is closer to the low polarity end of the polarity scale. The low polarity or less polar solvent can be a solvent selected from the group that includes acetone, ethanol, beta-phenethylamine, 2-ethoxyethanol, dioxane, methyl ethyl ketone, methyl n-propyl ketone, methyl acetate, methyl isobutyl ketone, chloroform, tetrahydrofuran, n-propanol, methyl isoamyl ketone, ethyl acetate, 2-methoxyethylacetate, isobutyl alcohol, n-butyl acetate, 2-butanol, 2-propanol, 1-butanol, ethylene dichloride, dichloromethane, ethyl ether, o-dichlorobenzene, chlorobenzene, benzene, o-xylene, m-xylene, p-xylene, methyl tertiary-butyl ether, toluene, carbon tetrachloride, trichloroethylene, n-butyl chloride, hexadecane, nonane, cyclohexane, trimethylpentane, petroleum ether, iso-hexanes, hexane, heptane, cyclopentane, trichlorotrifluoroethane, and pentane. In this example, the solvent for the detector reagent is preferably the least polar solvent in which the detector reagent has adequate solubility. The solvent for the detector reagent can be selected from the following list of solvents acetic acid, water, aqueous buffer solution with a pH in the range 2–12, dimethylsulfoxide, N-methylpyrrolidone, N,N-dimethyl acetamide, N,N-dimethyl formamide, propylene carbonate, acetonitrile, 2-methoxyethanol, diethylcarbonate, pyridine, methanol, acetone, ethanol, beta-phenethylamine, 2-ethoxyethanol, dioxane, methyl ethyl ketone, methyl n-propyl ketone, methyl acetate, methyl isobutyl ketone, chloroform, tetrahydrofuran, n-propanol, methyl isoamyl ketone, ethyl acetate, 2-methoxyethylacetate, isobutyl alcohol, n-butyl acetate, 2-butanol, 2-propanol, 1-butanol, ethylene dichloride, dichloromethane, ethyl ether, o-dichlorobenzene, chlorobenzene, benzene, o-xylene, m-xylene, p-xylene, methyl tertiary-butyl ether, toluene, carbon tetrachloride, trichloroethylene, n-butyl chloride, hexadecane, nonane, cyclohexane, trimethylpentane, petroleum ether, iso-hexanes, hexane, heptane, cyclopentane, trichlorotrifluoroethane, and pentane.

If the solvent in which the detector reagent is dissolved is too polar, then when the detector reagent solution is applied to the analyte spot on the TLC plate, the analyte spot may tend to become enlarged and the resulting analyte density in the spot may be decreased. Thus, for polar sorbent materials, where the solvent for the detector reagent solution is too polar, the detection sensitivity for the analyte of interest may be reduced.

For polar media such as ion-exchange TLC plates, in general, the solvent for the detector reagent should be an aqueous solution in which the predominant component is water.

For non-polar adsorbents such as reversed phase TLC plates, in general, the solvent or solvents included in a solvent mixture for dissolving the analytes as well as the detector reagent are preferably selected from solvents close to the high polarity end of the polarity scale. Solvents that can be selected for preparing a highly polar solvent mixture for use with reversed phase TLC plates include water,

methanol, N,N-dimethylformamide, acetonitrile, acetic acid, acetone, pyridine, ethanol, dioxane, chloroform, isopropanol, ethyl acetate, tetrahydrofuran, and n-propanol.

It is also possible to conduct a micro spot test by first applying the detector reagent to the TLC plate and allowing the solvent for the detector reagent to evaporate. Then, the solution containing the analyte is applied to the sorbent media or TLC plate. If this procedure is used, a detector reagent should be selected that is insoluble (or has very low solubility) in the solvent that is used to dissolve the analyte. If the detector reagent has some solubility in the solvent for the analyte, in a positive test result, a ring of the indicator may form instead of a small spot. Consequently, the detection sensitivity for the analyte may be poorer. Attractive advantages of applying the detector reagent to a TLC plate prior to applying the solutions containing the analyte(s) include (a) the liquid detector reagent solutions do not need to be prepared just prior to the test, (b) the required detector reagents can be pre-deposited at different locations on the same TLC media prior to on-site testing, and (c) the actual on-site testing steps are reduced to the microcapillary deposition of solutions containing the analytes and visual observation of the results.

To insure a long shelf life-stability for pre-deposited reagents, the prepared plates are preferably stored in dry condition, in an atmosphere that is free of moisture, and in which air has been displaced by an inert gas such as nitrogen or argon.

Most of the analytes that are separated and then detected using thin layer chromatography (TLC) or paper chromatography (PC), including, but not limited to, those detected using either a chromogenic or fluorogenic visualization reagent (often referred to as TLC or PC spray reagent), should be capable of being detected with high detection sensitivity using the micro spot test methodology described herein.

An illustrative methodology of carrying out the methods of the present invention is provided below:

A sample suspected of containing the analyte methylphosphonic acid is prepared by forming an acetone eluate from a polyester wipe. A microcapillary tube is used to draw up about 1 microliter of a solution containing the analyte and the end of the capillary tube is touched to a piece of a chromatographic sorbent medium, such as a TLC plate or medium. The analyte solution wets the sorbent layer by capillary action. Afterward, the TLC medium was allowed to dry and a drop of Bromcresol Green reagent was added. This caused a small yellow spot to be produced within a background of a large dark blue spot. This indicates that the acid has been retained near the spotting point due to its strong interaction with the chromatographic sorbent material. Since the analyte collects in a small area near the spotting point as a result of capillary action and being adsorbed into the TLC media, it is possible to detect minute quantities of the analyte. When the technique is used of bringing a microcapillary tube into contact with the surface of a TLC media, the analyte solution will exit from the microcapillary tube by capillary action. If the microcapillary tube is not kept in contact with the surface of the chromatographic media, a droplet larger than the diameter of the microcapillary tube may form. When a droplet larger than the diameter of the microcapillary tube forms and then comes in contact with the thin layer chromatographic media, the solution will wet a larger area and the analyte will not be as concentrated in a compact spot. Consequently, the detection sensitivity of the test may be poorer.

While this micro spot methodology has particular application to the detection of chemical warfare agents and examples of such applicability are given to demonstrate such applicability, the methodology is likewise applicable for use in conjunction with or in place of other thin layer chromatographic tests for environmental pollutants, contaminants, and hazards. Table 2 contains a list of compounds that are representative of the Priority 1 Analytes that can be detected with the processes of the present invention. This list represents a number of analytes that might be expected to be found during an on-site chemical weapons verification inspection. It will be understood by those of ordinary skill in the art that those analytes not specifically mentioned but known are also included herein and that the analyses for these analytes would be handled by the same methodology as analytes that are listed.

TABLE 2

PRIORITY 1 ANALYTES	
COMPOUND	SYNONYM
ethyl N,N-dimethylphosphoramidocyanate	GA
Isopropyl methylphosphonofluoridate	GB
Pinacolyl methylphosphonofluoridate	GD
Cyclohexyl methylphosphonofluoridate	GF
O-ethyl S-(2-diisopropylamino)ethyl methylphosphonothiolate	VX
bis(2-chloroethyl)sulfide	HD
bis[2-(2-chloroethylthio)ethyl]ether	T
2-chlorovinylchloroarsine	L
Methylphosphonic difluoride	DF
ethyl 2-(diisopropylamino)ethyl methylphosphonite	QL
Isopropyl methylphosphonic acid	IMPA
Pinacolyl methylphosphonic acid	PMPA
Cyclohexyl methylphosphonic acid	CMPA
Methylphosphonofluoric acid	MPFA
Methylphosphonic dichloride	DC
S-(2-diisopropylamino)ethyl methylphosphonothioic acid	EA 2192
ethyl methylphosphonic acid	EMPA
O-ethyl methylphosphonothioic acid	EMPTA
1,4-dithiane	DITHIANE
2-chlorovinylarsenious oxide	L-OXIDE
Methylphosphonic acid	MPA

According to another aspect of the invention, there is provided a method whereby a more specific indication of the analytes can be achieved by using two or more micro spot tests in combination. By using a series of spot tests, the user is able to accumulate evidence for or against the presence of a Priority 1 Analyte in the sample without actually identifying any of the specific chemical components of the sample. This is important because the acceptance of the on-site screening procedures by the chemical industry may ultimately depend on methodologies that minimize or eliminate the need for unnecessarily subjecting chemical samples to sophisticated, and potentially more intrusive, analytical methods.

If a sample unknown gives positive tests for one or more Priority 1 Analyte, TLCs could be used to determine if the suspect sample is a mixture, and to obtain R<sub>f</sub> value(s) of the suspect analyte(s). For example, the TLC can be used to show the relative positions (from which R<sub>f</sub> values are obtained) for spots resulting from, for example, phosphonic acids and dithiane. The data can be obtained using a procedure similar to that developed by Sass and Ludemann for the separation of phosphonic acids, see J. of Chromatography, 187, 447-452 (1980), the contents of which are incorporated herein by reference. It is also noteworthy to mention that the shape of a spot on the TLC media and the rate at which the spot becomes colored when contacted by the visualizing

reagent may also help to indicate which analyte is present. For example, a characteristic of the EMPTA spot is that it produces spots that have a long tail. Another characteristic of the EMPTA spot is that it changes color, going rapidly from colorless to brown when the TLC media is exposed to iodine vapor. While the micro spot test data may not be sufficient to identify the components of the unknown sample (which nonetheless is a desirable feature for screening tests), it is clear that the methods of the present invention can provide a considerable amount of evidence for the presence (or absence) of Priority 1 Analytes in a suspect sample.

Table 3 contains data that exemplify how three of the micro spot tests can be used in combination to accumulate presumptive evidence for the presence of several different Priority 1 Analytes. The sample unknown for the micro spot tests is one that would contain one of the following Priority 1 Analytes: MPA, EMPA, IMPA, PMPA, EMPTA and dithiane. For example, the data in Table 3 indicates that the response patterns from the three different spot tests can be used to distinguish dithiane and EMPTA from each other, and from MPA and several alkyloxy methylphosphonic acids that are also Priority 1 Analytes. A positive test result with the Bromcresol Green Test indicates that an acidic analyte, which could be MPA, or alkyloxy methylphosphonic acid, is in the sample. If the positive test with Bromcresol Green is combined with positive tests with TCNQ (7,7,8,8-tetracyanoquinodimethane) and gold chloride/NaOH, the response pattern could indicate that EMPTA may be present, but not dithiane, MPA or the alkyloxy methylphosphonic acids. A positive test with Bromcresol Green in combination with negative tests with TCNQ and gold chloride/NaOH indicates that a sample might contain MPA or one or more alkyloxy methylphosphonic acids, but not EMPTA or dithiane. Negative tests with Bromcresol Green and TCNQ combined with a positive gold chloride/NaOH test indicate that the sample may contain dithiane, but none of the phosphorus acids. Detection specificity is further improved when two or more tests are used in combination because different tests for the same analyte have different interference profiles.

TABLE 3

Results of Micro Spot Tests for Some Priority 1 Analytes			
Analyte	Reagent(s) for Micro Spot Test		
	Bromcresol Green	TCNQ	Gold Chloride/NaOH
MPA	+	-	-
EMPA	+	-	-
IMPA	+	-	-
PMPA	+	-	-
EMPTA	+	+	+
DITHIANE	-	-	+

Without additional data from other tests, however, the three spot tests used to obtain the data for Table 3, will not indicate if the unknown is a single substance or a mixture, and they will not indicate which phosphorus-containing acids may be present in the sample. However other tests could be used to provide more definitive results.

The following non-limiting examples further serve to illustrate the invention.

## EXAMPLE 1

## Micro Spot Test for Methylphosphonic Acids Using Bromocresol Green

## Detection Principle

A positive test is the appearance of a yellow spot in a larger blue spot. The control (analyte level=0) and negative tests are indicated by a blue spot that does not contain a yellow center. The color change that is observed in a positive test is due to the difference in the pH of the analyte (pH 3 and above) and the solid support (pH>5). At pH 3.8 and below, bromocresol green is yellow, and at pH 5.4 and above it is blue. The method described herein could be used for detecting other organic acids as well.

DETECTOR REAGENT: Bromocresol Green (0.04% in ethanol).

## Procedure for Preparing the Detector Reagent

The bromocresol green reagent is available from Aldrich Chemical Company. Transfer 2 ml of the reagent into a 3-ml plastic dropping bottle, replace the tip and screw on the cover.

SOLVENT FOR THE ANALYTE: An organic solvent (e.g. acetone, dichloromethane, hexane)

PREFERRED SOLID SUPPORT: MK6F Silica Gel 60A Glass Backed TLC Sheets, Clifton, N.J.

## Analytes Detected with this Test

Methylphosphonic acid (MPA), methylphosphonofluoridic acid (MPFA), ethyl methylphosphonic acid (EMPA), isopropyl methylphosphonic acid (IMPA), pinacolyl methylphosphonic acid (PMPA), cyclohexyl methylphosphonic acid (CMPA), O-ethyl methylphosphonothioic acid (EMPTA). It should be apparent to those skilled in the art that other acids that collect in a small spot when applied in solution to a TLC support, such as organic acids, should be capable of being detected in a similar manner.

## Detection Limit for Chemical Weapons Convention Analytes

MPA, EMPA, IMPA, PMPA, CMPA, EMPTA, and MPFA are detectable at the 100 ng level (i.e. when a 1 microliter aliquot of an acetone solution containing 0.01%, or more of analyte is spotted on preferred solid support using a microcapillary tube).

EMPA, IMPA, CMPA, PMPA and EMPTA are also detectable at the 10 ng level (i.e. when a 1 microliter aliquot of a dichloromethane solution containing 0.001% or more of analyte is spotted on the preferred solid support using a microcap).

## Equipment and Materials

- Locking forceps or spotting bulb assembly for holding microcap (e.g. cat. # 20-99, Analtech, Newark, Del.)
- Microcap, 1-microliter (e.g. cat. # 20-01, Analtech)
- Dropping Bottle, 3-ml capacity, (e.g. cat. # 211630, Wheaton, Millville, N.J.)
- 0.04% Bromocresol Green in Ethanol, cat. # B-7382, Sigma Chemical Co., St. Louis, Mo.
- MK6F Silica Gel 60A Glass Backed TLC Sheets or equivalent (e.g. cat. # 4861-1 10, Whatman Inc., Clifton, N.J.)
- Acetone (e.g. cat. # GC60032-4, Baxter Healthcare Corp., Burdick and Jackson Div., Muskegon, Mich.)
- Pencil

## Procedure

- Score a 1x3 inch TLC plate into twelve 0.5x0.5 inch sections with a pencil.
- Lock a 1-microliter microcap in the tip of the locking forceps.
- Place tip of microcap in a sample of pure acetone (or other solvent for the test) and wait a few seconds for the solvent to be drawn by capillary action to fill the microcap.

- Place the tip of the microcap in contact with the silica gel surface of the solid support near the center of one of the 0.5x0.5 inch sections. This is the 'control' (analyte level=0) spot.

- Wait a few seconds for the solvent to evaporate.

- Using a new microcap, for each sample, spot a different sample solution in each of the remaining 0.5x0.5 inch sections of the plate and allow the solvent to evaporate.

- Using the dropping bottle, add 1 drop of the bromocresol green to each spot.

- Observe the plate for the appearance of positive tests. A positive test is indicated by the appearance of a small yellow spot in a large green (wet) or blue (dry) spot. A positive detection signal appears within 1-2 seconds and the colors remain stable for at least several hours.

## Purpose and Applications

This micro spot test method is suitable for use as a field test for detecting analytes containing a phosphonic acid group. It provides evidence for or against the presence of chemical weapons convention (or CWC) analytes in a sample.. The test can be used alone or in conjunction with other micro spot tests that detect other functional groups in the sample. When two or more micro spot tests are used in combination, the detection specificity for target chemical weapons convention analyses is increased compared with the result of a single test.

## EXAMPLE 2

## Micro Spot Test for O-Ethyl Methylphosphonothioic Acid (EMPTA Using 7,7,8,8-Tetracyanoquinodimethane (TCNQ)

## Detection Principle

A positive test is the appearance of a blue spot in a larger pale yellow spot. If fresh reagent is not used, however, the reagent spot may be green instead of pale yellow. (With high analyte levels, the center of the blue spot may be bleached so that a white spot appears instead of a blue spot). The control (analyte level=0) and negative tests are indicated by a pale yellow spot that does not contain a blue spot in the center. The color change that is observed in a positive test is due to a sulfhydryl group in the analyte converting the 7,7,8,8-tetracyanoquinodimethane (TCNQ) reagent into a highly colored free radical. In a positive test, the color change occurs within 1 or 2 seconds after applying the TCNQ reagent.

DETECTOR REAGENT: 7,7,8,8-tetracyanoquinodimethane (2.5% in acetone)

## Procedure for Preparing Detector Reagent

In a 3-ml plastic dropping bottle place 5 mg of TCNQ. Add 2 ml of acetone. Place the dropping bottle tip in place and screw on the cap. Swirl until all of the TCNQ reagent dissolves.

SOLVENT FOR THE ANALYTE: Acetone, dichloromethane, or hexane.

PREFERRED SOLID SUPPORT: MK6F Silica Gel 60A Glass Backed TLC Sheets, cat.# 4861-110, Whatman Inc.

## Analytes Detected with this Test

O-ethyl methylphosphonothioic acid (EMPTA) as well as other materials containing phosphonothioic acid groups, sulfhydryl groups and other TCNQ free radical precursors. Detection Limits for Chemical Weapons Convention Analytes

EMPTA is detectable at the 10 ng level (i.e. when a 1-microliter aliquot of an acetone solution containing 0.001% or more of analyte is spotted on preferred solid supports using a microcapillary tube).

**Equipment and Materials**

same as Example 1 except that 7,7,8,8-Tetracyanoquinodimethane (e.g. cat. # B-7382, Sigma Chemical Co., St. Louis, Mo.) was used instead of the Bromcresol Green.

**Procedure**

same as Example 1 except that 7,7,8,8-Tetracyanoquinodimethane was used in step 7. In this example, the observation step required observing the plate for the appearance of positive tests which was indicated by the appearance of a small blue spot in a large yellow spot. The reagent spot may be green if fresh reagent is not used. A positive detection signal appears within 1-2 seconds and the colors remain stable for at least several hours.

**Purpose and Applications**

same as in Example 1.

**EXAMPLE 3**

Micro Spot Test for O-Ethyl  
Methylphosphonothioic Acid (EMPTA), 1,4-Dithiane, Bis(2-chloroethyl)sulfide (HD), and Bis[2-(2-ethylthio)ethyl]ether (T) Using Gold Chloride and Sodium Hydroxide

**Detection Principle**

In this test, two detector reagents are added in sequence. The first reagent is an aqueous solution of gold chloride. The second reagent is an aqueous solution of sodium hydroxide. It is believed that the first reagent forms a brown complex with compounds containing a thioether, phosphonothioic acid group or a sulfhydryl group. The second reagent, aqueous sodium hydroxide, probably hydrolyzes the complex thereby forming gold hydroxide, which is unstable, and decomposes to gold oxide. A purplish black spot (gold oxide) in a yellow background signifies a positive test. This color change occurs at the location where the sample was spotted on the solid support. Small black speckles may also appear in the test spot. The small speckles, which occur randomly in the reagent spot should be ignored.

**Detector Reagents**

1. Aqueous 4% Gold Chloride Solution

2. Aqueous 2N Sodium Hydroxide

**Procedure for Preparing Detector Reagents****Reagent #1**

Place hydrogen tetrachloroaurate trihydrate (1 g) in a 25 ml volumetric flask and add water to the mark. Allow the solution to stand for 1 week. Place 2 ml of the solution in 2-ml plastic dropping bottle. Replace the plastic tip and screw the cover on tightly.

**Reagent #2**

Place sodium hydroxide (8.0 grams) in a 100 ml volumetric flask. Add approximately 75 ml of water and swirl until the sodium hydroxide dissolves. Allow the solution to cool to room temperature. Add water to the mark.

**SOLVENT FOR THE ANALYTE:** An organic solvent (e.g. acetone, dichloromethane, hexane)

**PREFERRED SOLID SUPPORT:** MK6F Silica Gel 60A Glass Backed TLC Sheets, cat. # 4861-110, Whatman Inc.

**Analytes Detected with this Test**

O-ethyl methylphosphonothioic acid (EMPTA), Bis(2-chloroethyl)sulfide (HD), Bis[2-(2-ethylthioethyl)]ether (T), 1,4-dithiane and other compounds containing a thioether, a phosphonothioic acid group or a sulfhydryl group.

**Detection Limits for Chemical Weapons Convention Analytes**

EMPTA is detectable at the 1 ng level (i.e. when a 1-microliter aliquot of an dichloromethane or hexane solu-

tion containing 0.0001% or more of analyte is spotted on the solid support using a microcap).

Dithiane is detectable at the 10 nanogram level; i.e. when a 1-microliter aliquot of a dichloromethane or hexane solution containing 0.0001% or more of the analyte is spotted on the solid support using a microcap.

Bis[2-(2-ethylthio)ethyl]ether (T) is detectable at the 10 ng level when it is applied to the thin-layer chromatographic media in hexane solution. Bis(2-chloroethyl)sulfide is detectable at the 100 ng level when it is applied to the thin-layer chromatographic media in dichloromethane solution.

**Equipment and Materials**

Same as those used in Example 1 except that hydrogen tetrachloroaurate (III) trihydrate (e.g. cat. # 24,459-7, Aldrich Chemical Co.) and sodium hydroxide (#22146-5, Aldrich Chemical Co.) were used instead of the Bromcresol Green.

**Procedure**

The same first six steps of example 1 were followed. Thereafter,

7. Using the dropping bottle, add 1 drop of the gold chloride reagent to each spot.

8. Wait two minutes. [EMPTA can be detected down to the 10 ng level at this point. Therefore it can be distinguished from the other analytes (that require base-step 9)].

9. Using the dropping bottle, add 1 drop of the sodium hydroxide solution to each spot.

10. Observe the plate for the appearance of a positive test. A positive test is indicated by the appearance of a small purplish black spot in a large pale yellow spot. A positive detection signal appears within 1-2 seconds and the colors remain stable for at least several hours. The test spot may contain dark speckles that appear randomly and with increasing frequency as the spot ages. These should be ignored.

**EXAMPLE 4**

Micro Spot Test for Bis(2-chloroethyl)sulfide (HD), Bis[2-(2-ethylthio)ethyl]ether (T), and Other Mustards (Including Nitrogen Mustards) Using 4-(4'-Nitrobenzyl)pyridine and Sodium Hydroxide

**Detection Principle**

In this test, two detector reagents are used in combination. The first reagent is a 2% solution of 4-(4'-nitrobenzyl)pyridine in an organic solvent such as denatured ethyl alcohol or toluene. The second reagent is an aqueous solution of sodium hydroxide. The thin-layer chromatographic media is heated after the 4-(4'-nitrobenzyl)pyridine is applied to the analyte spot. In the first reaction, heat accelerates the alkylation of 4-(4'-nitrobenzyl)pyridine by the analyte. Basification then results in a deprotonation reaction that produces a blue dye. A positive test response is a small dark blue or purple spot on a white or pale red background.

**Detector Reagents**

1. 4-(4'-Nitrobenzyl)pyridine (2%) in denatured ethanol (or toluene)

2. Aqueous 2N Sodium Hydroxide

**Procedure for Preparing Detector Reagents****Reagent #1**

Place 4-(4'-nitrobenzyl)pyridine (20 mg) in a 2-ml plastic dropping bottle. Add 1 ml of acetone (or toluene). Swirl until the solid dissolves. Replace the plastic tip and screw the cover on tightly.

**Reagent #2**

Place sodium hydroxide (8.0 grams) in a 100 ml volumetric flask. Add approximately 75 ml of water and swirl

until the sodium hydroxide dissolves. Allow the solution to cool to room temperature. Add water to the mark.

**SOLVENT FOR THE ANALYTE:** Acetone, dichloromethane, or hexane.

**PREFERRED SOLID SUPPORT:** MK6F Silica Gel 60A Glass Backed TLC Sheets, cat # 4861-110, Whatman, Inc. Analytes Detected with this Test

Bis(2-chloroethyl)sulfide (HD), Bis[2-(2-ethylthio)ethyl] ether (T), as well as other sulfur- and nitrogen-mustards.

Other alkylating agents (e.g. diethyl sulfate) will also be detected.

**Detection Limits for Chemical Weapons Convention Analytes**

Bis [2-ethylthio)ethyl] ether (T) is detectable at the 10 ng level when it is applied to the thin-layer chromatographic media in hexane solution. Bis(2-chloroethyl)sulfide is detectable at the 100 ng level when it is applied to the thin-layer chromatographic media in dichloromethane solution.

#### Equipment and Materials

Same as those used in Example 1 except that 4-(4-nitrobenzyl)pyridine (e.g. cat. # N1,420-4, Aldrich Chemical Co.) and sodium hydroxide (#22146-5, Aldrich Chemical Co.) instead of Bromcresol Green.

#### Procedure

The same first six steps of example 1 were followed. Thereafter,

7. Using a dropping bottle, add 1 drop of the 4-(4'-nitrobenzyl)pyridine reagent.

8. Place the thin-layer chromatographic plate on a hot plate set at 90 degrees Centigrade. Wait for two minutes.

9. Remove the thin-layer chromatographic plate from the hot plate and allow it to cool for 15-30 seconds.

10. Using the dropping bottle, add 1 drop of the sodium hydroxide solution to each spot.

### EXAMPLE 5

#### Micro Spot Test for 4-(Dimethyl)aminopyridine Using Dragendorff Reagent

##### Detection Principle

A positive test is the appearance of a red spot in a larger orange spot. The color change that is observed in a positive test is due to reaction of 4-(Dimethylamino)pyridine with acetic acid and potassium bismuth iodide (Dragendorff Reagent) to produce a red complex.

##### Detect or Reagent

Whatman Inc. Catalog No. 4911-107 Dragendorff TLC Visualization Reagent containing potassium iodide and bismuth subnitrate in acetic acid.

##### Procedure for Preparing the Detector Reagent

Add 2 ml of the detector reagent to a 3-ml dropping bottle. Place the dropping bottle tip in place and screw on the cap.

**SOLVENT FOR THE ANALYTE:** Acetone, dichloromethane, or hexane:

**PREFERRED SOLID SUPPORT:** MK6F Silica Gel 60A Glass Backed TLC Sheets, cat #4861-110, Whatman, Inc.

Analytes Detected with this Test

4-(Dimethylamino)pyridine.

##### Detection Limits for Analytes

4-(Dimethylamino)pyridine is detectable at the 10 ng level (i.e. when a 1 microliter aliquot of an acetone solution containing 0.001 % or more of analyte is spotted on the preferred solid support using a microcapillary tube.

##### Equipment and Materials

Same as Example 1 except that Dragendorff Reagent (e.g. cat. # 4911-107 from Whatman, Inc., Fairfield, N.J. 07004) was used instead of Bromcresol Green.

#### Procedure

Same as Example 1 except that Dragendorff Reagent was used in step 7. In this example, the observation step required observing the plate for the appearance of positive tests which was indicated by the appearance of a small red spot in a large orange spot. A positive detection signal appears within 1 minute and the colors remain stable at least for several hours.

#### Purpose and Applications

This micro spot test is suitable as a test for detecting analytes such as 4-(dimethylamino)pyridine, quaternary nitrogen compounds, and alkaloids.

As a further example, the invention has applicability for urine and drug testing and as a supplement to commercially available thin layer chromatographic (TLC) test kits. One such test kit has been marketed by Eastman Kodak Co. (Cat. No. 13125 Kodak Chromat/O/Screen 60 Analysis Kit for Alkaloids). With this kit, the tests are performed by spotting a TLC strip, eluting the strip with a solvent, and then detecting the substances present by spraying the TLC strip with a chromogenic reagent. The appearance of a spot of a particular color indicates the presence of an alkaloid. Using the Chromat/O/Screen TLC kit as an example, the kit could be designed to use the microspot tests of the present invention as a means for rapidly prescreening a large number of samples. The microspot tests would be performed simply by spotting the analyte solution several times using a microcap each time on a different TLC strip, and then adding a drop of chromogenic detector reagent to each spot. Different chromogenic reagents would be used for each class of chemicals that are of interest. Compared with conventional TLC methods, the microspot test methodology of the present invention would result in increased sample throughput and much higher sensitivity of detection. Throughput would be increased because the time consuming elution step conventional with TLC tests, would not be needed, and samples that did not provide any positive responses in the microspot tests would be quickly eliminated from further consideration. If necessary, it would still be possible to use more time consuming TLC methods, but there would be a need to test only those samples that gave a positive test.

### REFERENCES

#### Detector Reagents

1. Jungreis, E., Spot Test Analysis: Clinical, Environmental, Forensic, and Geochemical Applications, Chemical Analysis, Volume 141, Second Edition, John Wiley and Sons, Inc., NY, 1997.
2. Mohammad, A., and Tiwari, S., Thirty-Five Years of Thin-Layer Chromatography in the Analysis of Inorganic Anions, Separation Science and Technology, 30(19), 3577-3614 (1995).
3. Jork, Helmut, Editor, Thin-Layer Chromatography: Reagents and Detection Methods, Vol. 1b, Physical and Chemical Detection Methods, John Wiley & Sons, NY, 1994.
4. Green, F. J., The Sigma Aldrich Handbook of Stains, Dyes, and Indicators, Aldrich Chemical Co., 1990.
5. Novak, T. J., and Davis, P. M., Detection of Sulfur Mustards Using Spectrofluorometry, U.S. Pat. No. 5,032, 380, Jul. 16, 1991.
6. Jork, Helmut, Editor, Thin-Layer Chromatography: Reagents and Detection Methods, Vol. 1a, Physical and Chemical Detection Methods, John Wiley & Sons, NY, 1989.
7. Novak, T. J., 4,4'-Dithiodianil, U.S. Pat. No. 4,414,414, Nov. 8, 1983.

8. Sherma, J., Practice and Applications of Thin Layer Chromatography on Whatman KC<sub>8</sub>Reversed Phase Plates, TLC Technical Series, Volume 1 (1981), Whatman Inc., Clifton, N.J. 07014.
9. Sass, S., and Ludemann, W., J. of Chromatography, 187, 447-452 (1980).
10. Gasparic, J., and Churacek, J., Detection Reagents, Laboratory Handbook of Paper and Thin Layer Chromatography, pp 323-335, Ellis Horwood, Ltd, England, 1978.
11. Corporate Authors, E. Merck, Dyeing Reagents for Thin Layer and Paper Chromatography, E. Merck, Darmstadt, Federal Republic of Germany, 1975.
12. Corporate Authors, Eastman Kodak Company, Eastman TLC Visualization Reagents and Chromatographic Solvents, Kodak Publication No. JJ-5 (1973), Eastman Kodak Company, Rochester, N.Y. 14650.
13. Zweig, G., and Sherma, J., Editors, Detection Reagents for Paper and Thin-Layer Chromatography, CRC Handbook of Chromatography, Volume II, Section II.1, pp 103-189, CRC Press, 1972.
14. Sawicki, E., Engel, C. R., and Elbert, W. C., Talanta 14, 1169-1178 (1967).
15. Ruch, W. E., Editor, Chemical Detection of Gaseous Pollutants: An Annotated Bibliography, Ann Arbor Science Publishers, Inc., Ann Arbor, Mich., 1966.
16. Feigl, F., Spot Tests in Organic Analysis, 6th Edition, Elsevier Science, Ltd., 1966.
17. Bryant, F., Overell, B. T., Biochim et Biophys. Acta 10, 471-6 (1963).
18. Feigl, F., Spot Tests in Inorganic Analysis, 7th Edition, Elsevier Science, Ltd, 1958.
19. Epstein, J., Rosenthal, R. W., and Ess, R. J., Anal. Chem. 27, 1435-39 (1955).
20. Bregoff, H. M., Roberts, E., Delwiche, C. C., J. Biol. Chem. 205, 565 (1953).
21. Munier, R. Bull. Soc. Chim. Biol. 35, 1225 (1953).
22. Obermiller, M., Angew. Chem. 49, 162-164 (1936).
23. Witten, B., and Prostack, A., Sensitive Detector Crayons for Phosgene, Hydrogen Cyanide, and Lewisite, Anal. Chem. 29, 885-7 (1957).
24. Pheil, R. W., Crayon for the Detection of G-Agents, U.S. Pat. No. 2,929, 791, Mar. 22, 1960.
25. Sass, S., Ludemann, W. D., Witten, B., Fischer, V., Sisti, A. J., and Miller, J. I., Colorimetric Determination of Certain Organophosphorus Compounds and Acylating Agents—Use of Diisonitrosoacetone Reagent, Anal. Chem. 29, 1346-9 (1957).
26. Kramer, D. N., and Morin, R. D., Detection of G-Agents, U.S. Pat. No. 2,926,072, Feb. 23, 1960.
27. Gehauf, B., Epstein, J., Wilson, G. B., Witten, B., Sass, S., Bauer, V. E., Rueggeberg, W. H. C., Anal. Chem. 29, 276 (1957).
28. Gehauf, B., and Goldenson, J., Detection and Estimation of Nerve Gases by Fluorescence Reaction, Anal. Chem. 29, 276 (1957).
29. Brante, G., Iodine as a Means of Development in Paper Chromatography, Nature 163, 651-2 (1949).
30. Sokolowski, M., and Rozylo, J. K., TLC Analysis of Warfare Agents under Battlefield Conditions, Journal of Planar Chromatography 6, 467-71 (1993).
31. Munavalli, S., and Pannella, M., Thin-Layer Chromatography of Mustard and Its Metabolites, Journal of Chromatog., 437, 423-8 (1988).
32. Ellman, G. H., Arch. Biochem. Biophys., 82, 70-77 (1959).
33. Mikrochim Acta, 788 (1971); Ibid., 341 (1973).

34. Mikrochim Acta, 526 (1972).
  35. Chemical and Engineering News, pg 29, Aug. 1, 1994.
  36. Yoe, J. H., and Sarver, L. A., Organic Analytical Reagents, John Wiley, NY, 1941, pages 66-326.
  37. Maile, R. J., Fishesser, G. J., and Anderson, M. M., Thin-Layer Chromatography of Phosphonic Acid Derivatives, Journal of Chromatog. 132, 366-68 (1977).
  38. Reiner, M., ed., Standard Methods of Clinical Chemistry, Vol. 1, Academic Press, NY, 1953, page 84.
- What is claimed is:
1. A method of detecting the presence of an analyte wherein the analyte remains at the spot of application and is analyzed at this same spot, comprising the steps of:
    - placing the analyte in a solution where the solvent for the analyte consists of a non-aqueous solvent selected from the group of dimethylsulfoxide, N-methylpyrrolidone, N,N-dimethyl acetamide, N,N-dimethyl formamide, propylene carbonate, acetonitrile, 2-methoxyethanol, diethylcarbonate, pyridine, methanol, acetone, ethanol, beta-phenethylamine, 2-ethoxyethanol, dioxane, methyl ethyl ketone, methyl n-propyl ketone, methyl acetate, methyl isobutyl ketone, chloroform, tetrahydrofuran, n-propanol, methyl isoamyl ketone, ethyl acetate, 2-methoxyethylacetate, isobutyl alcohol, n-butyl acetate, 2-butanol, 2-propanol, 1-butanol, ethylene dichloride, dichloromethane, ethyl ether, o-dichlorobenzene, chlorobenzene, benzene, o-xylene, m-xylene, p-xylene, methyl tertiary-butyl ether, toluene, carbon tetrachloride, trichloroethylene, n-butyl chloride, hexadecane, nonane, cyclohexane, trimethylpentane, petroleum ether, iso-hexanes, hexane, heptane, cyclopentane, trichlorotrifluoroethane, and pentane;
    - placing the solution containing the analyte in a tube having an end portion with a microcapillary sized opening, so that when the tube is placed in contact with a chromatographic sheet having a surface layer formed of sorbent material selected from the group consisting of silica gel, high performance thin layer chromatography (HPTLC) silica gel, polysilicic acid, aluminum oxide, cellulose, polyamide, reversed phase silica Gel C<sub>2</sub> (dimethyl bonded), reversed phase silica gel C<sub>2</sub> (ethyl bonded), reversed phase silica gel C<sub>8</sub> (octyl bonded), reversed phase silica gel C<sub>18</sub> (octadecyl bonded), acetylated cellulose, silica gel modified with amino groups, silica gel modified with cyano groups, Kieselghur impregnated with hydrocarbons, anionic and cationic anion exchange resins, diethylaminoethyl cellulose, and mixtures of the listed sorbents, the solution containing the analyte is withdrawn from the end portion of the tube and onto the surface layer of the sorbent material by capillary action;
    - placing the end portion of the tube having the microcapillary sized opening in contact with the sorbent material so the solution is withdrawn from the tube by capillary action, the solvent being absorbed into the sorbent material and the analyte being separated from the solvent and adsorbed by the sorbent material at the place of contact of the end portion of the tube with the sorbent material;
    - placing a chromogenic detector reagent for the analyte on the sorbent material at the place of contact of the end portion of the tube with the sorbent material to detect the presence of the analyte separated from the solvent and concentrated in the sorbent material at the place of contact of the end portion of the tube with the sorbent material, whereby a chromogenic indicator is formed when the analyte is present in the sorbent material.



2. A method of screening a solution for an analyte that has been dissolved in a solvent to form the solution and for detecting the presence of the analyte when the solution is deposited in a surface layer of a sorbent material wherein the analyte remains at the spot of application and is analyzed at this same spot comprising the steps of:

placing the solution containing the analyte in a tube having an end portion forming a microcapillary sized opening in the end portion of the tube so that when the tube is placed in contact with the sorbent material, the solution containing the analyte in the tube is withdrawn from the end portion of the tube and into the sorbent material by capillary action;

placing the end portion of the tube forming the microcapillary sized opening in contact with the sorbent material so that the solution is withdrawn from the tube by capillary action, the solvent being absorbed into the sorbent material and the analyte being separated from the solvent and adsorbed by the sorbent material at the point of contact of the end portion of the tube with the sorbent material; and

placing a detector reagent for the analyte on the sorbent material at the place of contact of the end portion of the tube with the sorbent material to detect the presence of the analyte separated from the solvent and concentrated in the sorbent material at the place of contact of the end portion of the tube with the sorbent material.

3. The method of claim 2, wherein the diameter of the microcapillary sized opening has range of diameters of from about 0.05 to about 1.6 millimeters so that when the end portion of the tube is placed in contact with the sorbent material, the solution containing the analyte is withdrawn from the end portion of the tube by capillary action, separated from the solvent at the place where the end portion of the tube having the microcapillary opening contacts the sorbent material and adsorbed by the sorbent material.

4. The method of claim 2, wherein the sorbent material is formed of a polar material selected from the group consisting of silica gel, high performance thin layer chromatography (HPTLC) silica gel, polysilicic acid, and aluminum oxide and the solvent for the analyte is a non-aqueous solvent that is less polar than the sorbent material and selected from the group of acetone, ethanol, beta-phenethylamine, 2-ethoxyethanol, dioxane, methyl ethyl ketone, methyl n-propyl ketone, methyl acetate, methyl isobutyl ketone, chloroform, tetrahydrofuran, n-propanol, methyl isoamyl ketone, ethyl acetate, 2-methoxyethylacetate, isobutyl alcohol, n-butyl acetate, 2-butanol, 2-propanol, 1-butanol, ethylene dichloride, dichloromethane, ethyl ether, o-dichlorobenzene, chlorobenzene, benzene, o-xylene, m-xylene, p-xylene, methyl tertiary-butyl ether, toluene, carbon tetrachloride, trichloroethylene, n-butyl chloride, hexadecane, nonane, cyclohexane, trimethylpentane, petroleum ether, isohexanes, hexane, heptane, cyclopentane, trichlorotrifluoroethane, and pentane.

5. The method of claim 2, wherein the sorbent material comprises a thin layer chromatographic sheet having a surface layer of silica gel sorbent material and the solvent for the analyte is selected from the group consisting of acetone, dichloromethane, toluene, o-xylene, m-xylene, p-xylene, n-butyl chloride, cyclohexane, trimethylpentane, petroleum ether, heptane, cyclopentane, pentane and hexane.

6. The method of claim 2, further comprising the step of dissolving the detector reagent in a detector reagent solvent to form a detector reagent solution prior to the step of placing the detector reagent for the analyte on the sorbent

material at the place of contact of the end portion of the tube with the sorbent material and to indicate the presence of the analyte separated from the solvent and concentrated in the sorbent material at the place of contact of the end portion of the tube with the sorbent material.

7. The method of claim 6, wherein the sorbent material is a polar material selected from the group of silica gel, high performance thin layer chromatography (HPTLC) silica gel, polysilicic acid, aluminum oxide, cellulose, and polyamide, and the solvent for the analyte is selected from solvents having less polarity than the sorbent material and selected from the group consisting of acetone, ethanol, beta-phenethylamine, 2-ethoxyethanol, dioxane, methyl ethyl ketone, methyl n-propyl ketone, methyl acetate, methyl isobutyl ketone, chloroform, tetrahydrofuran, n-propanol, methyl isoamyl ketone, ethyl acetate, 2-methoxyethylacetate, isobutyl alcohol, n-butyl acetate, 2-butanol, 2-propanol, 1-butanol, ethylene dichloride, dichloromethane, ethyl ether, o-dichlorobenzene, chlorobenzene, benzene, o-xylene, m-xylene, p-xylene, methyl tertiary-butyl ether, toluene, carbon tetrachloride, trichloroethylene, n-butyl chloride, hexadecane, nonane, cyclohexane, trimethylpentane, petroleum ether, isohexanes, hexane, heptane, cyclopentane, trichlorotrifluoroethane, and pentane.

8. The method of claim 6, wherein the sorbent material is a chromatographic material selected from the group consisting of silica gel, high performance thin layer chromatography (HPTLC) silica gel, polysilicic acid, aluminum oxide, cellulose, polyamide, reversed phase silica Gel C<sub>2</sub> (dimethyl bonded), reversed phase silica gel C<sub>2</sub> (ethyl bonded), reversed phase silica gel C<sub>8</sub> (octyl bonded), reversed phase silica gel C<sub>18</sub> (octadecyl bonded), acetylated cellulose, silica gel modified with amino groups, silica gel modified with cyano groups, Kieselghur impregnated with hydrocarbons, anionic and cationic anion exchange resins, diethylaminoethyl cellulose, and mixtures of the listed sorbents and the solvent for the analyte is selected from the group consisting of acetic acid, water, aqueous buffer solution with a pH in the range 2-12, dimethylsulfoxide, N-methylpyrrolidone, N,N-dimethyl acetamide, N,N-dimethyl formamide, propylene carbonate, acetonitrile, 2-methoxyethanol, diethylcarbonate, pyridine, methanol, acetone, ethanol, beta-phenethylamine, 2-ethoxyethanol, dioxane, methyl ethyl ketone, methyl n-propyl ketone, methyl acetate, methyl isobutyl ketone, chloroform, tetrahydrofuran, n-propanol, methyl isoamyl ketone, ethyl acetate, 2-methoxyethylacetate, isobutyl alcohol, n-butyl acetate, 2-butanol, 2-propanol, 1-butanol, ethylene dichloride, dichloromethane, ethyl ether, o-dichlorobenzene, chlorobenzene, benzene, o-xylene, m-xylene, p-xylene, methyl tertiary-butyl ether, toluene, carbon tetrachloride, trichloroethylene, n-butyl chloride, hexadecane, nonane, cyclohexane, trimethylpentane, petroleum ether, isohexanes, hexane, heptane, cyclopentane, trichlorotrifluoroethane, and pentane.

9. The method of claim 6, wherein the detector reagent is selected from the group consisting of bromocresol green; 7,7,8,8-tetracyanoquinodimethane (TCNQ); gold chloride; gold chloride/NaOH solution; 4-(4'-nitrobenzyl)pyridine/NaOH; cholinesterase/indoxyl acetate; sodium pyrophosphate peroxide/aromatic amine; potassium bismuth iodide; 1,3-diisonitrosoacetone guanidinium salt; bis(diethylamino) benzophenone oxime; bis(diethylamino)benzophenone; bis (dimethylamino)thiobenzophenone; phenylazoformic acid 2-diphenylhydrazide; diphenylcarbazon; diphenylthiocarbazon; mercuric salt; diethyldithiocarbamic acid silver salt;

29

2,2'-dithiobis(5-nitropyridine); 5,5'-dithiobis(2-nitrobenzoic acid), Ellman's Reagent; molybdenum oxide in sulfuric acid; ammonium molybdate; iodine/starch; and sulfuric acid (4M); ammonium sulfate; ammonium cerium(IV)sulfate; ammonium iron(II)sulfate; cobalt(II)thiocyanate; palladium (II)chloride; potassium iodide plateate; sodium tetraphenyl boron; o-tolidine; and N,2,6-trichloro-p-benzoquinoneimine.

10. The method of claim 9, wherein the solvent for the detector reagent has less polarity than the sorbent material.

11. The method of claim 2, wherein sorbent material is formed of a chromatographic polar material and the solvent for the analyte is a non-aqueous solvent that has a lower polarity than the sorbent material.

12. The method of claim 2, further comprising the step of dissolving the detector reagent in a detector reagent solvent to form a detector reagent solution prior to the step of placing the detector reagent for the analyte on the sorbent material, wherein the sorbent material is a polar chromatographic material and the solvents for the analyte and the detector reagents have less polarity than the sorbent material.

13. The method of claim 2, wherein the sorbent material is formed of a non-polar material selected from the group of reversed phase silica Gel C<sub>2</sub> (dimethyl bonded), reversed phase silica gel C<sub>2</sub> (ethyl bonded), reversed phase silica gel C<sub>8</sub> (octyl bonded), reversed phase silica gel C<sub>18</sub> (octadecyl bonded), acetylated cellulose, silica gel modified with amino groups, silica gel modified with cyano groups, and Kieselghur impregnated with hydrocarbons and the solvent for the analyte is an aqueous solvent mixture containing solvents from the group comprising water, methanol, N,N-dimethylformamide, acetonitrile, acetic acid, acetone, pyridine, ethanol, dioxane, chloroform, isopropanol, ethyl acetate, tetrahydrofuran, and n-propanol.

14. The method of claim 2, wherein the sorbent material is formed of an ion-exchange material selected from the group of anion exchange resin, cation exchange resin and diethylaminoethylcellulose and the solvent for the analyte comprises water.

15. A method of detecting the presence of an analyte wherein the analyte remains at the spot of application and is analyzed at this same spot, comprising the steps of:

placing a detector reagent for the analyte on a chromatographic layer of a sorbent material selected from the group consisting of silica gel, high performance thin layer chromatography (HPTLC) silica gel, polysilicic acid, aluminum oxide, cellulose, polyamide, reversed phase silica Gel C<sub>2</sub> (dimethyl bonded), reversed phase silica gel C<sub>2</sub> (ethyl bonded), reversed phase silica gel C<sub>8</sub> (octyl bonded), reversed phase silica gel C<sub>18</sub> (octadecyl bonded), acetylated cellulose, silica gel modified with amino groups, silica gel modified with cyano groups, Kieselghur impregnated with hydrocarbons, anionic and cationic anion exchange resins, and diethylaminoethyl cellulose;

placing the analyte in a solution where the solvent for the analyte is selected from the group of acetic acid, water, aqueous buffer solution with a pH in the range 2-12, dimethylsulfoxide, N-methylpyrrolidone, N,N-dimethyl acetamide, N,N-dimethyl formamide, propylene carbonate, acetonitrile, 2-methoxyethanol, diethylcarbonate, pyridine, methanol, acetone, ethanol, beta-phenethylamine, 2-ethoxyethanol, dioxane, methyl ethyl ketone, methyl n-propyl ketone, methyl acetate, methyl isobutyl ketone, chloroform, tetrahydrofuran, n-propanol, methyl isoamyl ketone,

30

ethyl acetate, 2-methoxyethylacetate, isobutyl alcohol, n-butyl acetate, 2-butanol, 2-propanol, 1-butanol, ethylene dichloride, dichloromethane, ethyl ether, o-dichlorobenzene, chlorobenzene, benzene, o-xylene, m-xylene, p-xylene, methyl tertiary-butyl ether, toluene, carbon tetrachloride, trichloroethylene, n-butyl chloride, hexadecane, nonane, cyclohexane, trimethylpentane, petroleum ether, iso-hexanes, hexane, heptane, cyclopentane, trichlorotrifluoroethane, and pentane;

placing the solution containing the analyte in a tube having an end portion formed with a microcapillary sized opening, so that when the tube is placed in contact with the sorbent material, the solution containing the analyte is withdrawn from the end portion of the tube and into the sorbent material by capillary action;

placing the end portion of the tube having the microcapillary sized opening in contact with the sorbent material at the spot where the detector reagent has been deposited on the sorbent layer so the solution is withdrawn from the tube by capillary action with the solvent being absorbed into the sorbent material and the analyte being separated from the solvent and adsorbed into the sorbent material at the point of contact of the end portion of the tube with the sorbent material, whereby a chromogenic indicator is formed when the analyte is present in the sorbent material.

16. A method of detecting the presence of an analyte in a solvent solution where the solution is deposited on a selected sorbent material by capillary action to cause the analyte to remain concentrated at the spot of deposition on the sorbent material wherein the analyte is analyzed at the same spot comprising the steps of:

placing the analyte in a solution where the solvent for the analyte consists of a non-aqueous solvent selected from the group of dimethylsulfoxide, N-methylpyrrolidone, N,N-dimethyl acetamide, N,N-dimethyl formamide, propylene carbonate, acetonitrile, 2-methoxyethanol, diethylcarbonate, pyridine, methanol, acetone, ethanol, beta-phenethylamine, 2-ethoxyethanol, dioxane, methyl ethyl ketone, methyl n-propyl ketone, methyl acetate, methyl isobutyl ketone, chloroform, tetrahydrofuran, n-propanol, methyl isoamyl ketone, ethyl acetate, 2-methoxyethylacetate, isobutyl alcohol, n-butyl acetate, 2-butanol, 2-propanol, 1-butanol, ethylene dichloride, dichloromethane, ethyl ether, o-dichlorobenzene, chlorobenzene, benzene, o-xylene, m-xylene, p-xylene, methyl tertiary-butyl ether, toluene, carbon tetrachloride, trichloroethylene, n-butyl chloride, hexadecane, nonane, cyclohexane, trimethylpentane, petroleum ether, iso-hexanes, hexane, heptane, cyclopentane, trichlorotrifluoroethane, and pentane;

placing the solution containing the analyte in a microcapillary sized tube having a volume of from about 0.1 to about 30.0 microliters and having an end portion with a microcapillary sized opening having a diameter of from about 0.05 to about 1.6 millimeters so that when the tube is placed in contact with a chromatographic sheet having a surface layer formed of sorbent material selected from the group consisting of silica gel, high performance thin layer chromatography (HPTLC) silica gel, polysilicic acid, aluminum oxide, cellulose, polyamide, reversed phase silica Gel C<sub>2</sub> (dimethyl bonded), reversed phase silica gel C<sub>2</sub> (ethyl bonded), reversed phase silica gel C<sub>8</sub> (octyl bonded), reversed phase silica gel C<sub>18</sub> (octadecyl bonded), acetylated

31

cellulose, silica gel modified with amino groups, silica gel modified with cyano groups, Kieselghur impregnated with hydrocarbons, anionic and cationic anion exchange resins, diethylaminoethyl cellulose, and mixtures of the listed sorbents, the solution containing the analyte is withdrawn from the end portion of the tube and onto the surface layer of the sorbent material by capillary action;

placing the end portion of the tube having the microcapillary sized opening in contact with the sorbent material so the solution is withdrawn from the tube by capillary action, the solvent being absorbed into the sorbent material and the analyte being separated from the

32

solvent and the analyte being concentrated and fixed at the point contact of the end portion of the tube with the sorbent material; and

placing a chromogenic detector reagent for the analyte on the sorbent material at the point of contact of the end portion of the tube with the sorbent material to detect the presence of the analyte that is concentrated in the sorbent material at the point of contact of the end portion of the tube with the sorbent material, whereby a chromogenic indicator is formed when the analyte is present in the sorbent material.

\* \* \* \* \*